



Line 1

Line 2

Patient ID

Clinician _____ Style _____ Qty _____

Patient Weight _____ Shoe Info: Youth Mens Womens Shoe Size: _____ Shoe Width: Narrow Std Wide

SHELL MODIFICATIONS

| | | | | | | | |
|-----------------------|-----------------------|------------|-------|-----------------------------|-----------------------|--------|------|
| Shell Rigidity | Flexible | Semi-Rigid | Rigid | Shell Length | Met | Sulcus | Full |
| Heel Cutout | L R B/L | | | 1 st Ray Cut Out | L R B/L | 45 ° | 60 ° |
| Flexible Morton's Ext | L R B/L | | | Semi-Rigid Morton's Ext | L R B/L | | |
| | Thickness: Consistent | 2.5mm | | | Thickness: Consistent | 2.5mm | |
| Reverse Morton's Ext | L R B/L | | | Heel Seat | 1/8" 1/4" 3/8" | | |
| | Thickness: Consistent | 2.5mm | | | 1/2" 5/8" 3/4" | | |
| Arch Modification | + 1/8" - 1/8" | | | Medial Flange | L R B/L | | |
| | | | | | Low Med High | | |
| | | | | Latelal Flange | L R B/L | | |
| | | | | | Low Med High | | |

POSTING OPTIONS

| | | | | | |
|----------------------------|-----------------------|---------------|---|-------------------------------|------------|
| Heel Lift (0 - 8 mm) | L _____ mm | R _____ mm | Medial Heel Skive (2mm 4mm 6mm 8mm) | L _____ mm | R _____ mm |
| | | | Lateral Heel Skive (2mm 4mm 6mm 8mm) | L _____ mm | R _____ mm |
| Rearfoot Post (0° - 6°) | ____ Lab Specify | | Forefoot Post (0° - 6°) | ____ Balance to Rearfoot | |
| | L ____ ° Varus Valgus | | | ____ Extrinsic ____ Intrinsic | |
| | R ____ ° Varus Valgus | | | L ____ ° Varus Valgus | |
| | | | | R ____ ° Varus Valgus | |
| Rearfoot Post Undercut | Narrow 20 ° | Standard 15 ° | | | |
| | Wide 10 ° | | | | |

TOP COVERS

| | | | | | |
|--------------------|---------------------|--------------------|--------------------|--------------------------------|--------------------|
| Top Cover Material | EVA 35 | ____ 1/8 ____ 1/16 | Mid Layer Material | EVA 35 | ____ 1/8 ____ 1/16 |
| | Micro-suede | ____ 1/24 | | Poly-U | ____ 1/8 ____ 1/16 |
| | Neolon | ____ 1/8 ____ 1/16 | | Poly-U w/Scrim | ____ 1/8 ____ 1/16 |
| | Plastazote | ____ 1/8 | Top Cover Length | ____ Met ____ Sulcus ____ Full | |
| | Poly-U (w/ vinyl) | ____ 1/8 ____ 1/16 | | | |
| | Ucolite Perforated) | ____ 1/8 ____ 1/16 | Bottom Cover | ____ Yes ____ No | |
| | Vinyl Black | ____ 1/24 | | | |

ACCOMMODATIONS

| | | | |
|--------------------|-------------|-----------------|-------------|
| Heel Cushion | L R B/L | Met Bar | L R B/L |
| Morton's Ext (EVA) | L R B/L | Dancer's Pad | L R B/L |
| Heel Spur Pad | L R B/L | Met Pads (1/8") | L R B/L |
| Neuroma Pad | L 1 2 3 4 5 | Balance Pad | L 1 2 3 4 5 |
| | R 1 2 3 4 5 | | R 1 2 3 4 5 |

NOTES
